


Helium-filled Soap Bubble Generator

for large volume PIV & Particle
Tracking



LaVision now developed in cooperation with DLR Göttingen and TU-Delft a new seeding generator and nozzle design for 100 or more nozzles in parallel. The bubbles have a highly constant size of down to 300 μm and are 10000 times brighter than standard oil aerosol seeding particles. With a relaxation time of 11 μs they follow the flow perfectly, even at high velocities and flow gradients [2]. Large measurement volumes can be filled with a high seeding concentration, recorded and processed with LaVision's time-resolved **Tomographic PIV** or **Shake-the-Box** particle tracking systems.



Bubbles leaving nozzle at a rate of 40000/sec.

Micro-processor controlled
operation

The micro-processor inside the controller device automatically manages the air, helium and soap flow rates with start, stop and pause modes. All maintenance procedures such as cleaning and filling are carried out at a push button. A large soap reservoir enables continuous operation of up to 20 h. Remote control is possible with LAN connection to any remote device with standard internet browser.

Multi-nozzle arrays

The nozzles are separated from the controller at a distance of up to 20 m. An aerodynamically shaped Linear Nozzle Array (LNA) with up to 20 nozzles has been designed for e.g. wind-tunnel operation, typically placed inside the settling chamber. Several LNAs can be connected to a single controller for flexible configurations and seeding large measurement volumes of more than 1 m^3 .



LaVisionUK Ltd

2 Minton Place / Victoria Road
Bicester, Oxon / OX26 6QB / United Kingdom
E-Mail: sales@lvision.com / www.lvisionuk.com
Phone: +44-(0)-870-997-6532 / Fax: +44-(0)-870-762-6252

LaVision GmbH


Anna-Vandenhoeck-Ring 19
D-37081 Göttingen / Germany
E-Mail: info@lvision.com / www.lvision.com
Tel. +49-(0)551-9004-0 / Fax +49-(0)551-9004-100

LaVision Inc.


211 W. Michigan Ave. / Suite 100
Ypsilanti, MI 48197 / USA
E-mail: sales@lvisioninc.com / www.lvisioninc.com
Phone: (734) 485 - 0913 / Fax: (240) 465 - 4306

Example

Thermal air convection in 1 cubic meter



Volume: $800 \times 800 \times 1200 \text{ mm}^3$,
illuminated with LED spots



Particle tracks and vorticity contours
calculated with Shake-the-Box (STB),
courtesy D. Schanz, A. Schröder, F. Huhn
DLR Göttingen

Specifications

- ▶ bubble size: down to $300 \mu\text{m}$
- ▶ 10000 brighter compared to usual oil aerosol seeding
- ▶ LED illumination possible in many applications
- ▶ neutrally buoyant, relaxation time $\tau=11 \mu\text{s}$
- ▶ bubble lifetime of several minutes
- ▶ aerodynamically shaped Linear Nozzle Array (LNA) with up to 20 nozzles
- ▶ large measurement volume $> 1 \text{ m}^3$
- ▶ high seeding concentration: 1300 bubbles per cm^3 at nozzle exit
- ▶ remote control via LAN: start, stop, pause, etc.

References:

- [1] Bosbach J, Kühn M, Wagner C (2009) Large scale particle image velocimetry with helium-filled soap bubbles. *Exp Fluids* 46:539–547
- [2] Scarano F, Sina Ghaemi S, Carlo Alp Caridi G, Bosbach J, Dierksheide U, Sciacchitano A (2015), On the use of helium-filled soap bubbles for large-scale tomographic PIV in wind tunnel experiments, *Exp Fluids* 56:42

Data provided by LaVision are believed to be true. However, no responsibility is assumed for possible inaccuracies or omissions. All data are subject to change without notice.

Jun-17

LaVisionUK Ltd

2 Minton Place / Victoria Road
Bicester, Oxon / OX26 6QB / United Kingdom
E-Mail: sales@lvision.com / www.lvisionuk.com
Phone: +44-(0)-870-997-6532 / Fax: +44-(0)-870-762-6252

LaVision GmbH

Anna-Vandenhoeck-Ring 19
D-37081 Göttingen / Germany
E-Mail: info@lvision.com / www.lvision.com
Tel. +49-(0)551-9004-0 / Fax +49-(0)551-9004-100

LaVision Inc.

211 W. Michigan Ave. / Suite 100
Ypsilanti, MI 48197 / USA
E-mail: sales@lvisioninc.com / www.lvisioninc.com
Phone: (734) 485 - 0913 / Fax: (240) 465 - 4306